

SHAKERdeluxe™



The new Shaker Deluxe is a simple and effective mucus clearance device designed to exporate trapped mucus, helping prevent infections and strengthening pulmonary defences.

The new adjustable and re-movable mouth piece design improves the usability and patient experience. Simply exhaling through the device creates a vibration in the airways which travels deep into the lung and shakes mucus loose, adjusting the angle of the device varies the intensity of the vibrations. Suitable for patients with COPD, Cystic Fibrosis Emphysema, Chronic Bronchitis, Bronchiectasis and Asthma.

- Respiratory physiotherapy device for mobilising pulmonary secretions such as mucus and catarrh
- Uses the latest innovative design technology
- Produces vibrations in the chest cavity to 'shake' stubborn mucus loose
- Aids expectoration, increases lung efficiency, helping to reduce breathlessness and fatigue
- Adjustable mouthpiece gives more user flexibility i.e. use sitting, standing or lying down
- Affordable, easy to use single user model, ideal for home healthcare or on the go
- 100's of thousands of Shaker Products sold worldwide



The Shaker Deluxe can be reused safely by a single person and ideally should not be shared.

www.powerbreathe-sa.co.za

**NOW AVAILABLE IN
SOUTH AFRICA**

Shaker Deluxe is not a toy. This product is designed to be used for breathing exercises only. Any other use is not recommended. The material in the manual is for information purposes only. POWERbreathe International Ltd makes no representations or warranties with respect to the manual or with respect to the products described. No claims are made or implied in the use or results of equipment. POWERbreathe International Ltd shall not be liable for any damages, direct or incidental, consequential or special, arising out of the use of this material or the products described. Shaker Deluxe is a patented and design protected product. All rights reserved. Specifications may change without notice due to manufacturers continuous programme of development. The POWERbreathe logo type is a registered trademark of POWERbreathe Holdings Ltd. All POWERbreathe product names are trademarks or registered trademarks of POWERbreathe Holdings Ltd. All details are correct at time of going to press. E & OE Ref: 3361

For more information:
Telephone: 011 440 2119 / 083 635 9567 email: shelelela@telkomsa.net

**POWER[®]
breathe**

SHAKERclassic™



The Shaker Classic is a simple and effective mucus clearance device designed to exporate trapped mucus, helping prevent infections and strengthening pulmonary defences.

Simply exhaling through the device creates a vibration in the airways which travels deep into the lung and shakes mucus loose, adjusting the angle of the device varies the intensity of the vibrations. Suitable for patients with COPD, Cystic Fibrosis Emphysema, Chronic Bronchitis, Bronchiectasis and Asthma.

- Offers simple and convenient relief from excessive mucus
- Can be used for chronic conditions such as COPD, asthma, emphysema and acute problems like chesty coughs, flu and bronchitis
- Adjustable mouthpiece gives more user flexibility i.e. use sitting, standing or lying down
- Affordable, easy to use single user model which can be used at home or on the move
- 100's of thousands of Shaker Products sold worldwide



The Shaker Classic can be reused safely by a single person and ideally should not be shared.

www.powerbreathe-sa.co.za

**NOW AVAILABLE IN
SOUTH AFRICA**

Shaker Classic is not a toy. This product is designed to be used for breathing exercises only. Any other use is not recommended. The material in the manual is for information purposes only. POWERbreathe International Ltd makes no representations or warranties with respect to this manual or with respect to the products described. No claims are made or implied in the use or results of equipment. POWERbreathe International Ltd shall not be liable for any damages, direct or incidental, consequential or special, arising out of the use of this material or the products described. Shaker Classic is a patented and design protected product. All rights reserved. Specifications may change without notice due to manufacturers continuous programme of development. The POWERbreathe logo type is a registered trademark of POWERbreathe Holdings Ltd. All POWERbreathe product names are trademarks or registered trademarks of POWERbreathe Holdings Ltd. All details are correct at time of going to press. E & OE Ref: 3361

For more information:
Telephone: 011 440 2119 / 083 635 9567 email: shelelela@telkomsa.net

**POWER
breathe**